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1. Yaziji H, Janckila AJ, Lear SC, Martin AW, Yam LT 1995 Immunohistochemical detection of tartrate-resistant acid phosphatase in non-hematopoietic human tissues. *Am J Clin Pathol* 104:397-402.
2. Lam KW, Li CY, Yam LT, Desnick RJ 1981 Comparison of the tartrate-resistant acid phosphatase in Gaucher's disease and leukemic reticuloendotheliosis. *Clin Biochem* 14:177-181.
3. Janckila AJ, Neustadt DH, Nakasato YR, Halleen JM, Hentunen T, Yam LT 2002 Serum tartrate-resistant acid phosphatase isoforms in rheumatoid arthritis. *Clin Chim Acta*. 320:49-58.
4. Halleen JM, Alatalo SL, Janckila AJ, Woitge HW, Seibel MJ, Väänänen HK 2001 Serum Tartrate-resistant acid phosphatase is a specific and sensitive marker of bone resorption. *Clin Chem* 47:597-600.
5. Hannon RA, Clowes JA, Eagleton AC, Al Hadari AA, Eastell R, Blumsohn A 2004 Clinical performance of immunoreactive tartrate resistant acid phosphatase isoform 5b as a marker of bone resorption. *Bone* 34:187-194.
6. Halleen JM, Alatalo SL, Suominen H, Cheng S, Janckila AJ, Väänänen HK 2000 Tartrate-resistant acid phosphatase 5b: a novel serum marker of bone resorption. *J Bone Miner Res* 15:1337-1345.
7. Janckila AJ, Takahashi K, Sun SZ, Yam LT 2001 Tartrate-resistant acid phosphatase isoform 5b as serum marker for osteoclastic activity. *Clin Chem* 47:74-80.
8. Albertini B, Casez JP, Cheneval JP, Wauters JP, Jaeger PH 2002 Low levels of 25-hydroxy-vitamin D3 (25OHD) and high serum tartrate resistant acid phosphatase 5b (TRAP 5b) activity in hemodialysis patients: relationship to secondary hyperparathyroidism. *J Am Soc Nephrol* 13:571A-572A.
9. Halleen JM 2003 Tartrate-resistant acid phosphatase 5B is a specific and sensitive marker of bone resorption (Review). *Anticancer Res* 23(2A):1027-1029.
10. Chu P, Chao TY, Lin YF, Janckila AJ, Yam LT 2003 Correlation between histomorphometric parameters of bone resorption and serum type 5b tartrate-resistant acid phosphatase in uremic patients on maintenance hemodialysis. *Am J Kidney Dis* 41:1052-1059.

11. Reichel H, Esser A, Roth HJ, Schmidt-Gayk H 2003 Influence of PTH Assay Methodology on Differential Diagnosis of Renal Bone Disease. *Nephrol Dial Transplant* 18:759-768.
12. Mehlhorn AT, Rechl H, Stemberger AW, Gradinger R 2002 Immunological and histochemical determination of TRAP 5b to measure osteoclast activity under antiosteolytic treatment. *Calcif Tissue Int* 70, Abstract P21.
13. Stepan JJ, Burckhardt P 2002 Serum activity of type 5b ACP and biochemical markers of type I collagen degradation in osteoporotic men with Klinefelter's syndrome treated with an intravenous ibandronate. *Calcif Tissue Int* 70, Abstract P79.
14. Halleen JM, Nenonen AM, Alatalo SL, Ivaska KK, Cheng S, Schmidt-Gayk H, Uusi-Rasi K, Heinonen A, Kannus P, Sievänen H, Väänänen HK 2003 Comparison of serum tartrate-resistant acid phosphatase 5b with other markers of bone turnover in monitoring alendronate therapy. *J Bone Miner Res* 18(Suppl.1), Abstract SA104.
15. Rissanen JP, Hentunen TA, Halleen JM 2003 Development and characterization of a novel human in vitro bone resorption assay useful for preclinical testing of drug candidates. *J Bone Miner Res* 18(suppl.1), Abstract SA256.
16. Terpos T, Viniou N, de la Fuente J, Meletis J, Voskaridou E, Karkantaris C, Vaiopoulos G, Palermos J, Yataganas X, Goldman J, Rahemtulla 2003 Pamidronate is superior to ibandronate in decreasing bone resorption, interleukin-6 and beta 2-microglobulin in multiple myeloma. *Eur J Haematol* 70:34-42.
17. Truniger R, Popp AWE, Perrelet R, Noesberger A, Lippuner K 2003 Effects of a 12 months risedronate treatment on biochemical markers of bone turnover and on serum osteoprotegerin and soluble receptor activator of NF-kB ligand in postmenopausal osteoporotic women. *J Bone Miner Res* 18(suppl. 1), Abstract SU341.
18. Voskaridou E, Terpos E, Spina G, Palmeros J, Rahemtulla A, Loutradi A, Loukopoulos D 2003 Pamidronate is an effective treatment for osteoporosis in patients with beta-thalassaemia. *Br J Haematol* 1234:730-737.
19. Alatalo SL, Halleen JM, Hentunen TA, Mönkkönen J, Väänänen HK 2000 Rapid screening method for osteoclast differentiation in vitro that measures tartrate-resistant acid phosphatase 5b activity secreted into the culture medium. *Clin Chem* 46:1751-1754.

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生产者名称：英国 Immunodiagnostic Systems Limited(IDS Ltd)  
生产者/生产场所地址：10 Didcot Way,Boldon Business Park,Boldon,Tyne&Wear,NE35 9PD,UK  
联系方式：+44 (0) 191 519 0660  
传真：+44 (0) 191 519 0760  
网址：www.idsltd.com

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