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1. Burgeson RE. New collagens, new concepts. Annu Rev Cell Biol 4:551-77;1988
2. Bonde M, Garnero P, Fledelius C, Qvist P, Delmas PD, Christiansen C. Measurement of bone degradation products in serum using antibodies reactive with an isomerized form of an 8 amino acid sequence of the C-telopeptide of Type I collagen. J Bone Miner Res 12:1028-1034; 1997
3. Garnero P, Gineyts E, Riou JP, Delmas PD. Assessment of bone resorption with a new marker of collagen degradation in patients with metabolic bone disease. J Clin Endocrinol & Metab 79:780-5; 1994
4. Bonde M, Qvist P, Fledelius C, Riis BJ, Christiansen C. Applications of an enzyme immunoassay for a new marker of bone resorption (CrossLaps<sup>®</sup>) - follow up on hormone

- replacement therapy and osteoporosis risk assessment. *J Clin Endocrinol & Metab* 80:864-8;1995
5. Pedersen BJ, Ravn P, Bonde M. Type I Collagen C-telopeptide Degradation products as Bone Resorption Markers. *J Clin LigandAssay* 21; 2: 118-127; 1998
  6. Rosenquist C, Fledelius C, Christgau S, Pedersen BJ, Bonde M, Qvist P, Christiansen C. SerumCrossLaps® ELISA. First application of monoclonal antibodies for measurement in serum of bone-related degradation products from C-terminal telopeptides of type I collagen. *Clin. Chem.*44;11: 2281-2289; 1998
  7. Christgau S, Rosenquist C, Alexandersen P, Bjarnason NH, Ravn P, Fledelius C, Herling C, Qvist P, Christiansen C. Clinical evaluation of the Serum CrossLaps® ELISA, a new assay measuring the serum concentration of bone-derived degradation products of type I collagen Ctelopeptides. *Clin Chem* 44:11; 2290-2300; 1998
  8. Bjarnason,N.H. & Christiansen,C. Early Response in Biochemical Markers Predicts Longterm Response in Bone Mass During Hormone Replacement Therapy in Early Postmenopausal Women. *Elsevier* 26, 561-569 (2000).
  9. Buclin T,Rochat MC,Burckhardt P,Azria M,Attinger M. Bioavailability and biological efficacy of a new oral formulation of salmon calcitonin in healthy volunteers. *J Bone Miner Res* 17, 8-1478 (2002).
  10. Christgau,S.; Bitsch-Jensen,O.; Bjarnason,N.H.; Gamwell Henriksen,E.; Qvist,P.; Alexandersen,P.; Bang Henriksen,D. Serum CrossLaps for Monitoring the Response in Individuals Undergoing Antiresorptive Therapy. *Elsevier* 26, 505-511 (2000).
  11. Delmas PD. Markers of bone turnover for monitoring treatment of osteoporosis with antiresorptive drugs. *Osteoporos Int* 11, S66-S76 (2000).
  12. Delmas,P., Hardy,P. & Dain,M.P. Monitoring Individual Response to Hormone Replacement Therapy with Bone Markers. *Elsevier* 26, 553-560 (2000).
  13. Fink,E.; Cormier,C.; Steinmetz,P.; Kindermans,C.; Bouc,Le Y.; Souberbielle,J.C. Differences in the Capacity of Several Biochemical Bone Markers to Assess High Bone Turnover in Early Menopause and Response to Alendronate Therapy. *Osteoporos. Int.* 11, 295-303 (2000).
  14. Ravn P,Neugebauer G,Christiansen C. Association between pharmacokinetics of oral ibandronate and clinical response in bone mass and bone turnover in women with postmenopausal osteoporosis. *Bone* 30, 320-324 (2002).
  15. Ravn,P. & Clemmesen,B. Biochemical Markers Can Predict the Response in Bone Mass During Alendronate Treatment in Early Postmenopausal Women. *Bone* 24, 237-244 (1997).
  16. Reginster JY,Henrotin Y,Christiansen C,Gamwell-Henriksen E,Bruyere,Collette J,Christgau S. Bone resorption in post-menopausal women with normal and low BMD assessed with biochemical markers specific for telopeptide derived degradation products of collagen type I. *Calcif Tissue Int* 69, 130-137 (2001).
  17. Rosen,H.N.; Moses,A.C.; Garber,J.; Iloputaife,I.D.; Ross,D.S.; Lee,S.L.; Greenspan,S.L. Serum CTX: A New Marker of Bone Resorption That Shows Treatment Effect More Often Than Other Markers Because of Low coefficient of Variability and Large Changes with Bisphosphonate Therapy. *Calcify Tissue Into* 66, 100-103 (2000).
  18. 18. S. Christgau, O. Bitch-Jensen, N. Hanover Bjarnason, E. Gamwell Henriksen, P. Qvist, P. Alexandersen and D. Bang Henriksen. Serum CrossLaps for Monitoring the Response in Individuals Undergoing Antiresorptive Therapy. *Bone* Vol. 26. No. 5 505-511 (2000)
  19. 19. Christiansen C, Tanko LB, Warming L, Moelgaard A, Christgau S, Qvist P, Baumann M, Wiczorek L, and Hoyle N. 2003. Dose dependent effects on bone resorption and formation of intermittently administered intravenous Ibandronate. *Osteoporos Int*.

20. 20. Chailurkit,L.O., S.Aunphongpuwanart, B.Ongphiphadhanakul, W.Jongjaroenprasert, S.Sae-Tung, and R.Rajatanavin. 2004. Efficacy of intermittent low dose alendronate in Thai postmenopausal osteoporosis. *Endocr Res* 30:29-36
21. 21. Qvist,P., M.Munk, N.Hoyle, and C.Christiansen. 2004. Serum and plasma fragments of C-telopeptides of type I collagen (CTX) are stable during storage at low temperatures for 3 years. *Clin Chim Acta* 350:167-173
22. 22. Reginster,J.Y., K.M.Wilson, E.Dumont, B.Bonvoisin, and J.Barrett. 2005. Monthly Oral Ibandronate is Well Tolerated and Efficacious in Postmenopausal Women: Results From the Monthly Oral Pilot Study. *J Clin Endocrinol Metab*
23. 23. Rosen,H.N., R.A.Parker, S.L.Greenspan, I.D.Iloputaife, L.Bookman, D.Chapin, I.Perlmutter, B.Kessel, P.Qvist, and M.Rosenblatt. 2004. Evaluation of Ability of Biochemical Markers of Bone Turnover to Predict a Response to Increased Doses of HRT. *Calcif Tissue Int*
24. 24. Tanko,L.B., Y.Z.Bagger, P.Alexandersen, J.P.Devogelaer, J.Y.Reginster, R.Chick, M.Olson, H.Benmammar, L.Mindeholm, M.Azria, and C.Christiansen. 2004. Safety and Efficacy of a Novel Salmon Calcitonin (sCT) Technology-Based Oral Formulation in Healthy Postmenopausal Women: Acute and 3-Month Effects on Biomarkers of Bone Turnover. *J Bone Miner Res* 19:1531-1538
25. 25. Warming,L., P.Ravn, D.Spielman, P.Delmas, and C.Christiansen. 2004. Trimegestone in a low-dose, continuous-combined hormone therapy regimen prevents bone loss in osteopenic postmenopausal women. *Menopause* 11:337-342.
26. 26. Warming,L., P.Ravn, T.Nielsen, and C.Christiansen. 2004. Safety and efficacy of drospirenone used in a continuous combination with 17beta-estradiol for prevention of postmenopausal osteoporosis. *Cimacteric* 7:203-111
27. 27. Warming,L., P.Ravn, and C.Christiansen. 2005. Levonorgestrel and 17beta-estradiol given transdermally for the prevention of postmenopausal osteoporosis. *Maturitas* 50:78-85.

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