

## Data Sheet

# HUMAN BONE SIALOPROTEIN (BSP)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 1110.1 / A 1110.2
<b>Immunogen:</b>	Recombinant human Bone Sialoprotein
<b>Synonyms:</b>	Bone sialoprotein 2, BSP2, Cell-binding sialoprotein, Integrin-binding sialoprotein
<b>Swiss-Prot No:</b>	P21815
<b>Gene Information:</b>	Gene name: IBSP GeneID: 3381
<b>Host:</b>	Goat
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Bone Sialoprotein
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA, RIA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	30 June 2016

### For research use only

**Publishing research using A 1110? Please let us know so that we can cite your publication as a reference.**

