

TECHNICAL SHEET

PRODUCT NAME: Anti-Giardia mAb trophozoite protein (clone G18)

Code: MT-16G18

Clone: G18

Isotype: IgG₁

Physical State: Liquid

Host Animal (source): Mouse

Specificity: G18 is specific for *Giardia lamblia* protein antigen.

Cross Reactivity: No cross reactivity with *Campylobacter spp., Clostridium difficile,*

Cryptosporidium parvum, Escherichia coli, Helicobacter pylori, Listeria monocytogenes, Salmonella spp., Shigella spp.,

Staphylococcus aureus, Entamoeba hystolitica.

Purification: Anti-Giardia mAb trophozoite protein (clone G18) has been

purified in a multi-step procedure.

≥95% pure by SDS-PAGE.

Storage Conditions: Storage at 2-8°C.

Health & Safety Information: This product contains sodium azide as preservative. Although the

amount is 0.09%, this substance is poisonous and hazardous and

it should be handled by trained staff only.

Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance

with the respective national regulations.

Presentation: PBS pH: 7.4, 0.09% Sodium Azide (NaN₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 11/01/2019



www.certest.es