

Data Sheet

HUMAN CHOLECYSTOKININ-RECEPTOR TYPE A [CCK-AR] (aa 180-195)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1129.1 / AE1129.2
Immunogen:	Synthetic human CCK-R Type A (aa 180-195) polyLysine conjugated (SNLVPFTKNNNQTANM)
Synonyms:	Cholecystokinin-1 receptor (CCK1-R)
Swiss-Prot No:	P32238
Gene Information:	Gene Name: CCKAR, CCKRA GeneID: 886
Host:	Rabbit
Matrix:	Serum
Specificity:	Human CCK-R Type A (aa 180-195)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	12 January 2018

For research use only

Publishing research using AE1129? Please let us know so that we can cite your publication as a reference.

