





## **Data Sheet**

## RAT GUANYLATE CYCLASE 2G (aa 1045-1057)

**ANTIBODY, POLYCLONAL** 

**Catalog no.:** AE1125.1 / AE1125.2

Immunogen: Synthetic rat Guanylate Cyclase 2G (aa 1045-1057) polyLysine-

conjugated (EFTEEEAKVPEIL)

**Synonyms:** Guanylyl cyclase receptor G (GC-G), Kinase-like domain-containing

soluble guanylyl cyclase (ksGC)

Swiss-Prot No: P55205

Gene Information: Gene Name: Gucy2g, Ksgc

GeneID:

Host: Rabbit
Matrix: Serum

**Specificity:** Synthetic rat Guanylate Cyclase 2G (aa 1045-1057)

Contents: 20 µl / 100 µl (lyophilized)

Resuspend in 20 µl / 100 µl aqua bidest.

**Known applications:** Western Blot (1:1000-2000)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive

controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)

Repeated thawing and freezing must be avoided

Last updated on: 17 June 2021

For research use only

Publishing research using AE1125? Please let us know so that we can cite your publication as a reference.

