

Data Sheet

HUMAN β -DEFENSIN 5 (aa 17-51)

ANTIBODY, POLYCLONAL

Catalog no.:	AE 1040.1 / AE 1040.2
Immunogen:	Synthetic human β -Defensin 5 (aa 17-51) (ESCKLGRGKCRKECLENEKPDGNCRLNFLCCRQRI)
Synonyms:	Beta-defensin 105, DEFB-5
Swiss-Prot No:	Q8NG35
Gene Information:	Gene Name: DEFB105A, BD5, DEFB105, DEFB5 and DEFB105B GenelD: 245908 and 504180
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human β -Defensin 5 (aa 17-51)
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	ELISA (1:2000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	19 August 2016

For research use only

Publishing research using AE1040? Please let us know so that we can cite your publication as a reference.

