

## Data Sheet

# HEMOGLOBIN

## ANTIBODY, POLYCLONAL

|                     |  |
|---------------------|--|
| <b>Catalog no.:</b> | A 7817.1 / A 7817.2  |
| <b>Immunogen:</b>   | Human Hemoglobin   |
| <b>Host:</b>        | Rabbit   |
| <b>Matrix:</b>      | Serum  |
| <b>Specificity:</b> | Human Hemoglobin   |
| <b>Contents:</b>    | 20 µl / 100 µl<br>Resuspend in 20 µl / 100 µl aqua bidest. |

**Known applications:** ELISA

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

**Store at:** 2-8 °C (lyophilized); - 20 °C (dissolved)  
Repeated thawing and freezing must be avoided

**Last updated on:** 19 August 2016

### For research use only

**Publishing research using A 7817? Please let us know so that we can cite your publication as a reference.**

