

Data Sheet

HUMAN GLUCAGON-LIKE PEPTIDE-2 [GLP-2]

ANTIBODY, POLYCLONAL

Catalog no.: A 6102.1 / A 6102.2

Immunogen: Synthetic human GLP-2, KLH-conjugated

(HADGSFSDEMNTILDNLAARDFINWLIQTKITD)

Synonyms: GLP-2, Glucagon (aa 126–158)

Swiss-Prot No: P01275

Gene Information: Gene Name: GCG

GenelD: 2641

Host: Rabbit **Matrix:** Serum

Specificity: Human GLP-2

Contents: 20 μl / 100 μl (lyophilized)

Resuspend in 20 µl / 100 µl aqua bidest.

Known applications: ELISA (1:3000)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should

determine optimal dilutions in their system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)

Repeated thawing and freezing must be avoided

Last updated on: 29 September 2015

For research use only

Publishing research using A 6102? Please let us know so that we can cite your publication as a reference.

