

Data Sheet

HUMAN C-PEPTIDE

ANTIBODY, POLYCLONAL, AFFINITY PURIFIED

Catalog no.: A 5512.1 / A 5512.2

Immunogen: C-Peptide conjugate EAEDLQVGQVKKKC-KLH

Synonyms: Insulin
Swiss-Prot No: P01308

Gene Information: Gene Name: INS

GeneID: 3630

Host: Sheep

Preparation: Affinity purified from sheep serum against target peptide and cross-

adsorbed on an affinity column of immobilised mouse immunoglobulin

Matrix: PBS pH7.4, 0.1% NaN₃

Specificity Human C-Peptide

Contents: 10μg / 100 μg (lyophilized)

Resuspend in 10 µl / 100 µl aqua bidest.

Known applications: ELISA

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)

Repeated thawing and freezing must be avoided

Last updated on: 17 January 2018

For research use only

Publishing research using A 5512? Please let us know so that we can cite your publication as a reference.

