



Data Sheet

HUMAN SECRETORY IMMUNOGLOBULIN A (sIgA)

ANTIBODY, POLYCLONAL

Catalog no.:	A 1083.1 / A 1083.2
Immunogen:	Human sIgA
Host:	Sheep
Matrix:	Serum
Specificity:	Human sIgA
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:5000) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	05 July 2016

For research use only

Publishing research using A 1083? Please let us know so that we can cite your publication as a reference.

