

Data Sheet

PORCINE DIAMINE OXIDASE (DAO)

ANTIBODY, POLYCLONAL

Catalog no.: A 1026.1 / A 1026.2

Immunogen: Porcine Diamine oxidase

Synonyms: Amiloride-sensitive amine oxidase, Amiloride-binding protein (ABP),

Histaminase, PK-DAO

Swiss-Prot No: Q9TRC7

Gene Information: Gene Name: AOC1, ABP1

GeneID: 100517436

Host: Rabbit **Matrix:** Serum

Specificity: Porcine DAO

Contents: 20 μl / 100 μl (lyophilized)

Resuspend in 20 µl / 100 µl aqua bidest.

Known applications: ELISA, Western Blot (1:5000)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should

determine optimal dilutions in their system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)

Repeated thawing and freezing must be avoided

Last updated on: 29 August 2017

For research use only

Publishing research using A 1026? Please let us know so that we can cite your publication as a reference.

