

TECHNICAL SHEET

PRODUCT NAME: Enterovirus EV71 recombinant protein

Code: MT-25ETV71, MT-25ETV71s1

Physical State: Liquid

Source: Escherichia coli

Description: Enterovirus EV71 recombinant VP1 protein expressed and

purified from Escherichia coli.

Appearance: Clear

Purity: >85%

Storage Conditions: Short term storage at 2-8°C. Long term storage at

-20°C. Avoid multiple/freeze thaw cycles by storing multiple

aliquots at -20°C.

Health & Safety Information: This product contains sodium azide as preservative. Although

the amount is 0.09%, this substance is poisonous and hazardous

and it should be handled by trained staff only.

Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in

accordance with the respective national regulations.

Presentation: Tris buffer with additives and salts, pH 8.0, 0.09% Sodium Azide

(NaN₃).

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 21/01/2019



<u>certest@certest.es</u> <u>www.certest.es</u>