

TECHNICAL SHEET

PRODUCT NAME: SC2	2 recombinant Spike Glycoprotein (S) UK variant
Code:	MT-25C19S.b.1.1.7, MT-25C19S.b.1.1.7s1
Physical State:	Liquid
Source:	Hek293 cells
Description:	SC2 recombinant Spike Glycoprotein (S) UK variant expressed and purified from <i>Hek293 cells</i> .
Appearance:	Clear
Purity:	>90%
Storage Conditions:	Short term storage at 2-8°C. Long term storage at -20°C. Avoid multiple/freeze thaw cycles by storing multiple aliquots at -20°C.
Shipping:	With cool packs.
Health & Safety Informat	tion: Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.
Presentation:	Tris buffer with additives and salts, pH: 8.00.

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 28/04/2021





Pol. Industrial Río Gállego II, Calle J, № 1, 50840, San Mateo de Gállego, Zaragoza (SPAIN) Tel. +34 976520354 <u>certest@certest.es</u> <u>www.certest.es</u>