

TECHNICAL SHEET

PRODUCT NAME: SARS-CoV-2 recombinant Nucleoprotein (NP) BR-P.1

Code: MT-25C19NP.p.1, MT-25C19NP.p.1s1

Physical State: Liquid

Source: Escherichia coli

Description: SARS-CoV-2 recombinant Nucleoprotein (NP) BR-P.1 expressed

and purified from Escherichia coli.

Appearance: Clear

Purity: >90%

Storage Conditions: Short term storage at 2-8°C. Long term storage at -20°C. Avoid

multiple/freeze thaw cycles by storing multiple aliquots at

-20ºC.

Shipping: With cool packs.

Health & Safety Information: Good Laboratory Practices should be followed when handling

this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in

BIOTEC

accordance with the respective national regulations.

Presentation: Tris buffer with additives and salts, pH: 8.00.

THIS PRODUCT IS INTENDED FOR RESEARCH USE ONLY.

Date: 28/04/2021

